



Plastic Tub Rail

Item #	Product Size	Color	UOM
12039	12"	White	1/bx

- Attractive, ergonomically-designed tub rail provides a comfortable, safe grip.
- Handle is perpendicular to tub.
- 13" in height.
- Tool-free, easy to install.
- Width adjusts from 2.75" to 6.25".
- Retail packaged.

specifications

Adjustable Clamp	2.75" (smallest) 6.25" (widest)
Height	13"
Weight	6 lbs.
Carton Shipping Weight	9 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Plastic Tub Rail

Item #	Product Size	Color	UOM
12044	16½"	White	1/bx

- Attractive, ergonomically-designed tub rail provides a comfortable, safe grip.
- Handle is perpendicular to tub.
- 16" in height.
- Tool-free, easy to install.
- Width adjusts from 3" - 6.25".
- Retail packaged.

specifications

Adjustable Clamp	3" (smallest) 6.25" (widest)
Height	16"
Carton Shipping Weight	7 lbs.
Weight	5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime