



**LIQUICELL**  
COMFORT EQUIPPED

FULL SIZE SCOOTER

Seat Color  
Black Vinyl

Body Colors  
Red  
Blue

# Ride

## GOLDEN COMPANION

4 wheel

The Companion four-wheel scooter offers the most leg room and foot room in its class with a modern, rugged design. Features include easy disassembly; an adjustable LED headlight; ergonomic control panel design with backlit battery gauge display; and wraparound Delta tiller. The Companion four-wheel is rugged, dependable and an outstanding value for a full-time, four-wheel scooter.



Rise. Roll. Rest. Ride.

**GOLDEN**  
technologies

YOUR LIFE IN MOTION

## Companion 4 wheel Specifications

<b>Model</b>	<b>GC440</b>
Medicare Code	K0807
Weight Capacity	350 lbs
Drive Wheels	Rear
Maximum Speed	4.5 mph
Operating Range <sup>1</sup>	15 miles
Ground Clearance to center deck	4.5"
Ground Clearance to transaxle/motor	2.5"
Turning Radius	53.5"
Type Batteries	2-U1
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	90 Amp Dynamic R Series
<b>Dimensions<sup>3</sup>:</b>	
Length	47.5"
Width	24"
Height (ground to top of back of seat)	36.5" - 39.5"
Ground to Top of Seat	23" - 26"
Ground to Top of Deck	6.25"
Deck to Top of Seat	17" - 20"
Number of Seat Height Positions	4
Size of Increments	1"
Front Axle to Rear Axle	35.5"
<b>Weight of Unit: (assembled)</b>	186 lbs
Front Half	63 lbs
Rear Half	48 lbs
Seat w/Arms	28 lbs
Batteries	22 lbs ea.
<b>Tires:</b>	Solid/Low Profile
Front	9"
Rear	10"
Rear Anti-Tippers	3"
<b>Standard Seat<sup>3</sup>:</b>	Stadium Style
Slide Seat	Standard
Back Height - no headrest	16"
Width x Depth	18" X 16"
Color	Black Vinyl
<b>Optional Seats<sup>3</sup>:</b>	Stadium Style
Back Height - no headrest	17.5"
Width	16", 18" or 20"
Depth	16" or 18"
Color	Charcoal Vinyl, Nubuck Vinyl
<b>Optional Seats<sup>2</sup> 3:</b>	Captain's Highback
Back Height - no headrest	21"
Back Height with headrest	29"
Width	18"
Depth	18"
Color	Coffee/Sand, Charcoal/Gray
<b>Seat Options:</b>	
Pneumatic Rise Seat Option	Yes
Electric Rise Seat Option	Yes
<b>Warranty:</b>	
Frame	Lifetime
Drive Train	2 Years
Electronics <sup>2</sup>	2 Years



Angle adjustable LED headlight creates a wide path of bright light



EZ Reach tiller adjustment




Patented infinite adjustable tiller



Quiet & maintenance free motor

## Standard Features on Every Model:

- Under carpet troubleshooting guide
- Delta tiller with auto lock mechanism
- EZ Reach tiller adjustment lever
- LiquiCell Comfort Equipped seat
- Seat slide adjustment lever
- Adjustable LED headlight



Blue Red



Disassembles easily!

<sup>1</sup> Battery range will vary due to rider weight, drive surface and terrain.

<sup>2</sup> Electronics warranty excludes batteries.

<sup>3</sup> There is a +1"/-1" tolerance on all dimensions.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.