

Mattress - Alternating Pressure Support Surface

Features

Inner spring :

- Cell-in-Cell design keeps internal bladders filled in case of a power failure, eliminating the need for a foam base.
- Removable, non-shear, anti-microbial fully enclosed cover encapsulates the entire 20 bladder mattress.
- CPR feature allows for rapid deflation of mattress system.
- Overall dimensions : 35" x 79" x 8".
- Weight : 12.4 lbs.
- Alternating cell depth : 8".

Powerloft PRO APM500 PUMP:

- High flow pump with 4-6 liters per minute.
- Soft/Firm Dial to obtain custom patient comfort pressure levels.
- Low pressure Alarm LED, an audio visual alarm indicating when the mattress pressure is very low (when mattress hose is disconnected or when pump stops working).
- Static function to enable cells from alternating.
- Hospital grade UL/CSA approved power cord.
- Low pressure Alarm mute function.
- Voltage Range : 110-120 / 220-240 VAC
- Frequency Range : 50 / 60 Hz
- Dimension : 12" x 7.5" x 4"
- Weight : 6.2 lbs