

## Deluxe 3-In-1 Steel commode

C311-4 Seat height: 16" - 21" Assembled 4/ctn

### Features

- Three-in-one commode reduces the number of items carried
- Durable plastic snap-on seat with lid
- Easily raised or lowered from 16" to 22" in 1" increments
- Plastic arms for added comfort
- Comes complete with commode bucket, cover and splash shield
- Nylon capped nut and bolt construction
- 300lb weight capacity



C311



C312

## Deep Seat 3-In-1 Steel Commode

C312-4 Seat height: 16" - 22" Assembled 4/ctn

### Features

- Deep seat provides extra comfort and support
- Three-in-one commode reduces the number of items carried
- Frame construction made by steel
- Plastic armrests for added comfort
- Comes with commode bucket, cover and splash shield

## Drop-Arm Steel Commode

C313-2 Seat height: 16.5" - 22.5" Assembled 2/ctn

### Features

- Easy to release arm mechanism allows users transfer themselves to and from commode
- Manufactured with 1" steel tube frame providing more strength and stability
- Powder coated frame is easy to maintain
- Padded armrests provide extra comfort and support
- Comes with commode pail with carry handle
- Limited lifetime warranty



C313



## Deluxe 3-In-1 Alum. Commode

C321-4 Seat height: 17" - 23" Assembled 4/ctn

### Features

- Three-in-one commode reduces the number of items carried
- Durable plastic snap-on seat with lid
- Easily raised or lowered from 17" to 23" in 1" increments
- Foam arms for added comfort
- Comes complete with commode bucket, cover and splash shield
- Nylon capped nut and bolt construction
- Frame completely assembled
- Aluminum construction
- Frame completely assembled
- 300lb weight capacity



C321



C322

## Deep Seat 3-In-1 Alum. Commode

C322-4, Seat height: 16.5" - 22.5" Assembled 4/ctn

### Features

- Deep seat provides extra comfort and support
- Three-in-one commode reduces the number of items carried
- Frame construction made by aluminum
- Foam arms for added comfort
- Comes with commode bucket, cover and splash shield
- Limited lifetime warranty

## Raised Toilet Seat with Lock

C401-1 1/bx

### Features

- Heavy Duty blow-molded plastic construction offers additional strength and durability
- Entire foam padded armrest provides more comfortable
- Improved locking device to ensure a safe and secure fit
- Limited lifetime warranty



C401

## Bariatric 3-In-1 Steel Commode

C314-2 Seat height: 16.5" - 22.5" 2/ctn

### Features

- Accommodates individuals up to 650lb
- Back is easy to install/ remove by wing nuts
- Comes with commode bucket and cover
- Extra durable frame construction by steel tubing
- PU armrests provide comfort and stability
- Limited lifetime warranty



C314

## Extra-Large Heavy Duty, Drop Arm Commode

C315-1, C316-1, C317-1 1/ctn

### Features

- 650lb (C315) / 850lb (C316) / 1000lb (C317) Weight Capacity
- Extra molded PU foam, extra wide and deep toilet seat
- All steel construction
- Easy-to-release drop arm mechanism allows for safe lateral transfer to and from commode
- PU armrests provide comfort and support
- Comes complete with commode bucket and cover
- Limited lifetime warranty



C315



C316



C317

Item	Width (Between Arms) (cm/in)	Outside Width (cm/in)	Seat Width (cm/in)	Seat Depth (cm/in)	Seat Height (cm/in)	Height (cm/in)	Weight (cm/in)	Weight Limit (kg/lb)	Weight Limit (kg/lb)
C314	71/28	78/31	35/14	43/17	42-57/16.5-22.5	67-82/26.5-32.5	16.1/35.5	<del>35.9/79</del>	300/650
C315	62/24.5	70/27.5	61/24	45/18	47/18.5	63/25	16.5/36.5	<del>20.7/45.5</del>	300/650
C316	62/24.5	70/27.5	61/24	53/21	47/18.5	63/25	16.8/37	<del>24.7/47</del>	385/850
C317	93/36.5	100/39.5	61/24	53/21	47/18.5	63/25	19.5/43	<del>24.5/54</del>	454/1000