



H101

## Trapeze

H101-1 1/ctn

### Features

- Mounts easily to metal frame headboards
- Chrome plated finish for two piece square steel tubes
- Easy to install without tools
- Easy to adjust height, length and horizontal position
- Adjustment levers on legs
- For use with homecare style beds

## Hydraulic Manual Patient Lift

H102-1 1/ctn

A Sling: Heavy-Gauge, Canvas, without helmet

B Sling: Deluxe Slings with Helmet

### Features

- Ideal for safe transfer of patients
- Easy to set up and disassemble
- High performance hydraulics raise/lower individuals
- Pump handle rotates from side to side for attendant convenience
- Padded swivel bar
- Invacare convertible type
- Limited lifetime warranty (frame)
- 1 Year limited warranty (pump)



H103-1 1/ctn

### Features

- Aluminum hydraulic pump with quick release valve
- Maximum 330lb (150kg) load capacity
- 2" round tube construction
- Power coating finish
- Castor: front 4" dia, rear 4" dia, with brake device



H103