



Overbed Table

T100-1 1/ctn

Features

- Knob adjustment for height from 26" - 38.5"
- Durable wood grain laminated top with T-mounted edge
- Swivel casters allow for easy maneuverability



T100

Overbed Table

T101-1 1/ctn

Features

- Deluxe overbed table provide the tabletop angle and height adjustments
- Knob adjustment for height from 26" - 38.5"
- Tabletop angle adjusts to 30 - 60 - 90 degrees in two directions
- Durable wood grain laminated top with T-mounted edge
- Swivel casters allow for easy maneuverability



T101

Auto-touch Overbed Table

T110-1 1/ctn

Features

- Light height adjustment touch on auto locking handle
- allows table to elevate freely
- Adjustment for height from 26" - 38.5"
- Durable wood grain laminated top with T-mounted edge
- Swivel casters allow for easy maneuverability



T110

Full Auto-touch Overbed Table

T111-1 Tabletop angle adjustable 1/ctn

Features

- Light height adjustment touch on auto locking handle allows table to elevate freely
- Adjustment for height from 26" - 38.5"
- Auto locking handle allows tabletop angle adjustment to 30 - 60 - 90 degrees
- Durable wood grain laminated top with T-mounted edge
- Swivel casters allow for easy maneuverability



T111