

Medicare Code

E0185

Features

- The Sierra™ Gel Overlay incorporates a series of gel bladders encapsulated between high-density polyurethane sculptured foam and support base.
- This combination creates a therapeutically advanced yet cost effective support surface.
- Optimal environment for prevention/treatment of pressure sores.
- Easy to clean, non-shear/low friction surface.
- Elastic corner straps to secure overlay to mattress.
- Base foundation is a high-density solid one-piece foam constructed core.
- Convoluted high-density TM, a soft nylon urethane coated OFX non-vapor permeable material.
- Bladders contain an aqueous non-toxic gel solution.
- Size hospital twin.
- Components conform to CAL TECH 117 and all infection control standards.
- 34" x 72" x 3"