



## Deluxe Folding Walker, One Button

W110-4 Adult, height: 32" - 39" 4/ctn

W111-4 Junior, height: 25" - 32" 4/ctn

### Features

- Single-button release allows quick, one-hand folding
- Crossbrace is located near hand grips to allow for full stride
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight
- 300lb weight capacity
- EZ fold for storage



W110 / W111



W120 / W121

## Deluxe Folding Walker, Two Button

W120-4 Adult, height: ~~32" - 39"~~ 4/ctn 33"-40"

W121-4 Junior, height: ~~25" - 32"~~ 4/ctn 26"-33"

### Features

- Easy to use red button release can be operated by fingers, palm, or side of hand
- Crossbrace is located near hand grips to allow for full stride
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight
- 300lb weight capacity
- EZ fold for storage

## Deluxe Folding Walker, Two Button with 5" wheels & glide tip

W124-4 Adult, height: ~~32" - 39"~~ 4/ctn 33"-40"

W125-4 Junior, height: ~~25" - 32"~~ 4/ctn 26"-33"

### Features

- Easy to use red button release can be operated by fingers, palm, or side of hand
- Crossbrace is located near hand grips to allow for full stride
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight
- 300lb weight capacity
- Includes 5" (Accessory# 010) front wheels and rear glide tip (Accessory# 011)
- EZ fold for storage



W124 / W125

## Deluxe Folding Walker, Universal size

- W130-4 Height adjustable from 28" - 38" One Button 4/ctn
- W132-4 Height adjustable from ~~28" - 38"~~ Two Button 4/ctn
- W133-4 Height adjustable from ~~28" - 38"~~ Two Button with 5" wheels

### Features

29"-39"

- Single-button release allows quick, one-hand folding
- Easy to use red button release can be operated by fingers, palm, or side of hand
- Crossbrace is located near hand grips to allow for full stride
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight
- 300lb weight capacity
- EZ fold for storage
- W133 is for W132+5" wheel foottube



## Reciprocal Walker, One Button

- W112 Height adjustable from ~~32" - 39"~~ 4/ctn

### Features

30.5"-36.5"

- Single-button release allows quick, one-hand folding
- Crossbrace is located near hand grips to allow for full stride
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight
- 300lb weight capacity
- EZ fold for storage
- Frame swivels enabling the frame to move with the user





## Bariatric Aluminum Folding Walker, Two Button

W140-2 Adult, height: 32" - 39" 2/ctn

W141-2 Junior, height: 25" - 32" 2/ctn



W140

### Features

- For Bariatric use
- Oversized frame design provides wider and deeper to accommodate individuals up to 500lb
- Precision design provides additional strength while adding minimal weight
- Limited lifetime warranty


## Bariatric Aluminum Folding Walker, Two Button with 5" wheels

W144-2 Adult, height: 32" - 39" 2/ctn

W145-2 Junior, height: 25" - 32" 2/ctn

### Features

- Same features as W140 / W141
- Includes 5" front wheels (Accessory# 010)
- Includes rear glide tips (Accessory# 011) allowing use on all surfaces
- Limited lifetime warranty



W144

## Side Walker (Hemi-Walker)

W150-2 height: ~~29.5" - 37"~~ 2/ctn

29.3"-36"

### Features

- Designed for individuals with the use of only one hand or arm
- Lighter than a walker and more stable than a cane
- Folds easily with one hand
- Ideal for users who do not require a walker but need a wider base for support
- Limited lifetime warranty



W150

## Accessories For Folding Walker



012 Walker basket



W144+Glide tip



W130+5" Wheel Foottube



011 Glide tip 1/pr



013 Glide tip 1/pr



010+5" Wheel Foottube 1/pr