

Lift Recliner Specifications



TRADITIONAL SERIES

| | Cambridge | Williamsburg | Pioneer |
|--|-----------|--------------|----------|
| Model Number | PR-401 | PR-747 | PR-643 |
| Size | Medium | Medium | Medium |
| Weight Capacity | 375 Lb. | 375 Lb. | 375 Lb. |
| Recline Positions | 3 | 3 | 3 |
| Standard Fabrics | 5 | 5 | 5 |
| Fabric Options | 28 | 32 | 32 |
| Chaise Seat | YES | YES | YES |
| Back Type | Biscuit | Button | Seam |
| Head Rest Cover | YES | YES | YES |
| Arm Covers | YES | YES | NO |
| Battery Back Up | YES | YES | YES |
| Overall Width | 33" | 32.5" | 29.5" |
| Overall Height | 41.5" | 44" | 41" |
| Width Between Arms | 20" | 21" | 21.5" |
| Floor to Top of Seat | 19" | 19.5" | 19" |
| Seat Depth | 20" | 21.5" | 21" |
| Seat to Top of Back | 28" | 28" | 28" |
| Seat to Top of Arm | 7" | 6" | 8" |
| Floor to Back Edge of Seat Raised | 28" | 28" | 27.5" |
| Floor to Front Edge of Seat Raised | 24" | 24" | 23" |
| Seat Angle Raised | 20 Deg. | 15 Deg. | 20 Deg. |
| Floor to Top of Foot Rest Reclined | 22" | 23" | 23" |
| Floor to Top of Back Reclined | 33" | 36" | 34" |
| Back Angle Closed | 110 Deg. | 100 Deg. | 115 Deg. |
| Back Angle Reclined | 145 Deg. | 130 Deg. | 140 Deg. |
| Foot Rest Angle Closed | 90 Deg. | 90 Deg. | 90 Deg. |
| Foot Rest Angle Opened | 175 Deg. | 180 Deg. | 180 Deg. |
| Distance Required from Wall Reclined | 20" | 18" | 18" |
| Extension Reclined from Back of Seat to end of Foot Rest | 41.5" | 42" | 40.5" |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism | YES | YES | YES |
| Warranty on Electrical Parts | 3 Yrs. | 3 Yrs. | 3 Yrs. |
| Pro Rated Years 4 Thru 7 | YES | YES | YES |
| Warranty on Mechanical Labor | 3 Yrs. | 3 Yrs. | 3 Yrs. |
| Warranty on Electrical Labor | 1 Yr. | 1 Yr. | 1 Yr. |
| Removable Back Standard | YES | YES | YES |
| Smart Tek System Standard | YES | YES | YES |
| Customizable | NO | NO | NO |
| Left Hand Control Option Available | YES | YES | NO |
| Extra Pocket Option Available | YES | YES | YES |
| Heat & Massage Option Available | YES | YES | YES |
| Foot Rest Extension Option Available | YES | YES | NO |
| Straight Lift Option Available | YES | YES | YES |
| User Height (suggested) | | | |



SIGNATURE SERIES

| | Winston | Regal | Royal | Space Saver |
|----------|----------|----------|----------|-------------|
| PR-410 | PR-751TY | PR-752 | PR-906 | |
| Large | Medium | Medium | Medium | |
| 375 Lb. | 375 Lb. | 375 Lb. | 250 Lb. | |
| 3 | 3 | 3 | 3 | |
| 4 | 4 | 4 | 4 | |
| 4 | 32 | 32 | 32 | |
| YES | YES | YES | YES | |
| Biscuit | 3 Pillow | 3 Pillow | Seam | |
| YES | YES | YES | YES | |
| NO | NO | YES | NO | |
| YES | YES | YES | YES | |
| 43" | 36" | 32.5" | 35.5" | |
| 46" | 44" | 44" | 43" | |
| 21.5" | 22" | 21" | 20" | |
| 21" | 20" | 19.5" | 22" | |
| 21" | 21.5" | 21" | 22" | |
| 29" | 28" | 28" | 27" | |
| 6.5" | 8" | 6" | 8" | |
| 29" | 28" | 27.5" | 30" | |
| 24.5" | 24" | 23.5" | 24" | |
| 20 Deg. | 18 Deg. | 15 Deg. | 24 Deg. | |
| 23.5" | 24" | 23.5" | 21.5" | |
| 36" | 33" | 33" | 37" | |
| 110 Deg. | 110 Deg. | 115 Deg. | 120 Deg. | |
| 158 Deg. | 140 Deg. | 145 Deg. | 140 Deg. | |
| 90 Deg. | 90 Deg. | 90 Deg. | 90 Deg. | |
| 180 Deg. | 180 Deg. | 180 Deg. | 180 Deg. | |
| 19" | 18" | 19" | 6" | |
| 45" | 42" | 42" | 44" | |
| YES | YES | YES | 1 Yr. | |
| 3 Yrs. | 3 Yrs. | 3 Yrs. | 1 Yr. | |
| YES | YES | YES | NO | |
| 3 Yrs. | 3 Yrs. | 3 Yrs. | 1 Yr. | |
| 1 Yr. | 1 Yr. | 1 Yr. | 1 Yr. | |
| YES | YES | YES | YES | |
| YES | YES | YES | YES | |
| NO | NO | NO | NO | |
| YES | NO | YES | YES | |
| YES | YES | YES | YES | |
| YES | YES | YES | YES | |
| NO | NO | YES | NO | |
| YES | YES | YES | NO | |



MAXICOMFORT

| | Relaxer | Comforter | Comforter | Comforter |
|-------------|-------------|-------------|-------------|-----------|
| PR-756MC | PR-505S | PR-505M | PR-505L | |
| Medium | Small | Medium | Large | |
| 375 Lb. | 300 Lb. | 375 Lb. | 375 Lb. | |
| MaxiComfort | MaxiComfort | MaxiComfort | MaxiComfort | |
| 4 | 4 | 4 | 4 | |
| 32 | 32 | 32 | 32 | |
| YES | YES | YES | YES | |
| 4 Pillow | Seam | Seam | Seam | |
| YES | YES | YES | YES | |
| NO | YES | YES | YES | |
| YES | YES | YES | YES | |
| 36.5" | 31.5" | 32.5" | 34.5" | |
| 44" | 40" | 42" | 44" | |
| 20" | 20" | 21" | 23" | |
| 19" | 19" | 19.5" | 20.5" | |
| 21" | 19.5" | 21" | 22" | |
| 29" | 25" | 27" | 28.5" | |
| 8" | 6" | 7" | 5.5" | |
| 26.5" | 27.5" | 27.5" | 29" | |
| 22" | 24" | 23.5" | 25" | |
| 15 Deg. | 15 Deg. | 16 Deg. | 16 Deg. | |
| 23" | 23" | 23" | 23" | |
| 21" | 32" | 32" | 33" | |
| 100 Deg. | 100 Deg. | 100 Deg. | 100 Deg. | |
| 174 Deg. | 174 Deg. | 174 Deg. | 174 Deg. | |
| 90 Deg. | 90 Deg. | 90 Deg. | 90 Deg. | |
| 180 Deg. | 180 Deg. | 180 Deg. | 180 Deg. | |
| 25" | 20" | 21" | 22" | |
| 43" | 39" | 42" | 43" | |
| YES | YES | YES | YES | |
| 3 Yrs. | 3 Yrs. | 3 Yrs. | 3 Yrs. | |
| YES | YES | YES | YES | |
| 3 Yrs. | 3 Yrs. | 3 Yrs. | 3 Yrs. | |
| 1 Yr. | 1 Yr. | 1 Yr. | 1 Yr. | |
| YES | YES | YES | YES | |
| YES | YES | YES | YES | |
| NO | NO | NO | NO | |
| YES | YES | YES | YES | |
| YES | YES | YES | YES | |
| YES | YES | YES | YES | |
| NO | NO | NO | NO | |
| NO | NO | NO | NO | |



COMFORTER SERIES

| | Comforter | Comforter | Comforter | Comforter | Comforter |
|------------|-----------|-----------|-----------|-----------|-----------|
| PR-501JP | PR-501S | PR-501M | PR-501L | PR-501T | |
| Jr. Petite | Small | Medium | Large | Tall | |
| 300 Lb. | 300 Lb. | 375 Lb. | 375 Lb. | 375 Lb. | |
| 3 | 3 | 3 | 3 | 3 | |
| 5 | 5 | 5 | 5 | 5 | |
| 32 | 32 | 32 | 32 | 32 | |
| YES | YES | YES | YES | YES | |
| Seam | Seam | Seam | Seam | Seam | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| 30.5" | 31.5" | 32.5" | 34.5" | 34.5" | |
| 39" | 40" | 42" | 44" | 47" | |
| 19" | 20" | 21" | 23" | 23" | |
| 18" | 19" | 19.5" | 20.5" | 22.5" | |
| 18.5" | 19.5" | 21" | 22" | 24" | |
| 22.5" | 25" | 27" | 28.5" | 30" | |
| 6.5" | 6" | 7" | 5.5" | 6" | |
| 26" | 27.5" | 27.5" | 29" | 31" | |
| 22" | 24" | 23.5" | 25" | 26" | |
| 16 Deg. | 15 Deg. | 16 Deg. | 16 Deg. | 18 Deg. | |
| 22" | 23" | 23" | 23" | 23" | |
| 32" | 32" | 32" | 33" | 36" | |
| 115 Deg. | 115 Deg. | 115 Deg. | 115 Deg. | 115 Deg. | |
| 145 Deg. | 145 Deg. | 145 Deg. | 145 Deg. | 145 Deg. | |
| 90 Deg. | 90 Deg. | 90 Deg. | 90 Deg. | 90 Deg. | |
| 180 Deg. | 180 Deg. | 180 Deg. | 180 Deg. | 180 Deg. | |
| 17" | 17" | 17" | 18" | 19" | |
| 38" | 39" | 42" | 43" | 45" | |
| YES | YES | YES | YES | YES | |
| 3 Yrs. | 3 Yrs. | 3 Yrs. | 3 Yrs. | 3 Yrs. | |
| YES | YES | YES | YES | YES | |
| 3 Yrs. | 3 Yrs. | 3 Yrs. | 3 Yrs. | 3 Yrs. | |
| 1 Yr. | 1 Yr. | 1 Yr. | 1 Yr. | 1 Yr. | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| NO | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| NO | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |



COMFORTER WIDES

| | Comforter | Comforter | Comforter | Comforter | Comforter |
|------------|-------------|-------------|-------------|-------------|-----------|
| PR-501S-23 | PR-501M-26D | PR-501L-26D | PR-501T-28D | PR-502 | |
| Wide | Medium | Large | Tall | Wide | |
| 375 Lb. | 500 Lb. | 500 Lb. | 500 Lb. | 700 Lb. | |
| 3 | 3 | 3 | 3 | Independent | |
| 4 | 4 | 4 | 4 | 4 | |
| 32 | 32 | 32 | 32 | 32 | |
| YES | YES | YES | YES | NO | |
| Seam | Seam | Seam | Seam | Seam | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| 34.5" | 38.5" | 38.5" | 41" | 45" | |
| 40" | 42" | 44" | 48" | 47" | |
| 23" | 27.5" | 27.5" | 29.5" | 33" | |
| 18" | 19.5" | 20.5" | 22.5" | 21" | |
| 19.5" | 21.5" | 22.5" | 24" | 22" | |
| 25" | 27" | 28.5" | 30" | 29" | |
| 6" | 7" | 5.5" | 6" | 7" | |
| 27.5" | 27.5" | 29" | 31" | 31" | |
| 24" | 23.5" | 25" | 26" | 27" | |
| 15 Deg. | 16 Deg. | 16 Deg. | 18 Deg. | 12 Deg. | |
| 23" | 23" | 23" | 23" | 21" | |
| 32" | 32" | 33" | 36" | 39" | |
| 115 Deg. | 115 Deg. | 115 Deg. | 115 Deg. | 110 Deg. | |
| 145 Deg. | 145 Deg. | 145 Deg. | 145 Deg. | 150 Deg. | |
| 90 Deg. | 90 Deg. | 90 Deg. | 90 Deg. | 90 Deg. | |
| 180 Deg. | 180 Deg. | 180 Deg. | 180 Deg. | 180 Deg. | |
| 17" | 17" | 18" | 19" | 20" | |
| 39" | 42" | 43" | 45" | 45" | |
| YES | YES | YES | YES | 2 Yrs. | |
| 3 Yrs. | 3 Yrs. | 3 Yrs. | 3 Yrs. | 1 Yr. | |
| YES | YES | YES | YES | NO | |
| 3 Yrs. | 3 Yrs. | 3 Yrs. | 3 Yrs. | 2 Yrs. | |
| 1 Yr. | 1 Yr. | 1 Yr. | 1 Yr. | 1 Yr. | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | NO | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | NO | |
| YES* | YES* | YES* | YES* | YES* | |



VALUE SERIES

| | Capri | Monarch | Monarch | Monarch Plus | Monarch Plus |
|----------|----------|----------|----------|--------------|--------------|
| PR-200 | PR-355M | PR-355L | PR-359M | PR-359L | |
| Medium | Medium | Large | Medium | Large | |
| 325 Lb. | 375 Lb. | 375 Lb. | 375 Lb. | 375 Lb. | |
| 2 | 3 | 3 | 3 | 3 | |
| 5 | 5 | 5 | 5 | 5 | |
| 5 | 32 | 32 | 32 | 32 | |
| NO | NO | NO | YES | YES | |
| Button | Seam | Seam | Seam | Seam | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| 31" | 32" | 34" | 32" | 34" | |
| 41" | 42" | 44" | 41.5" | 43" | |
| 20" | 21" | 23" | 20" | 23" | |
| 18" | 18.5" | 18.5" | 19.5" | 20.5" | |
| 19" | 21.5" | 22" | 21" | 22" | |
| 27" | 28" | 30" | 27.5" | 29.5" | |
| 7" | 7" | 7" | 6" | 5.5" | |
| 27" | 27" | 27" | 28" | 29" | |
| 22" | 23" | 23" | 24" | 25" | |
| 16 Deg. | 16 Deg. | 16 Deg. | 16 Deg. | 16 Deg. | |
| 23" | 23" | 22" | 23" | 23" | |
| 35" | 33" | 33" | 33" | 33" | |
| 110 Deg. | 110 Deg. | 110 Deg. | 110 Deg. | 110 Deg. | |
| 135 Deg. | 140 Deg. | 140 Deg. | 140 Deg. | 140 Deg. | |
| 90 Deg. | 90 Deg. | 90 Deg. | 90 Deg. | 90 Deg. | |
| 160 Deg. | 176 Deg. | 176 Deg. | 176 Deg. | 176 Deg. | |
| 11" | 17" | 18" | 17" | 18" | |
| 40" | 43" | 43" | 43" | 43" | |
| YES | YES | YES | YES | YES | |
| 1 Yr. | 2 Yrs. | 2 Yrs. | 2 Yrs. | 2 Yrs. | |
| NO | NO | NO | NO | NO | |
| 6 Mos. | 6 Mos. | 6 Mos. | 6 Mos. | 6 Mos. | |
| NO | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| NO | NO | NO | NO | NO | |
| YES | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |
| NO | YES | YES | YES | YES | |
| NO | YES | YES | YES | YES | |
| YES | YES | YES | YES | YES | |

