

Cirona™ 6300

Deep Vein Thrombosis Prevention Therapy System

Product Overview

The **Cirona™ 6300** is a portable DVT compression therapy pump that works by applying intermittent pressure to the patient's limbs affected by DVT. This therapy mimics the body's natural muscle contractions which can help increase blood circulations and decrease venous pressure. It may also reduce wound healing time and can aid in the treatment of stasis dermatitis, venous stasis ulcers, and chronic venous insufficiency. The **Cirona™ 6300** can be used in ambulatory, post-operative, and home care settings.

Features and Benefits:

- ✓ **Portable:** Suitable to use while on the go or at home
- ✓ **Convenient:** Comes with rechargeable batteries that may last for up to **8 hours** of treatment
- ✓ **Simple:** Easy to use one button operation
- ✓ **Comfortable:** Adjustable sleeves are made of soft, breathable fabric

Device Specifications:

Weight:	0.3kg / 0.66lbs
Dimensions:	65 x 130 x 30 (mm) 2.6 x 5.1 x 1.2 (in)
AC/DC Adapter:	AC Input: 100~240 Volts, AC 50/60Hz, 500mA MAX DC Output: 5V 2A (2 Plugs)
Operation Mode:	Continuous / 8 hours
Battery Type:	Lithium battery 3.7v 1350mAh (rechargeable)
Charging Time:	3 Hours
Default Pressure:	50 mmHg ± 5 mmHg
Rest Period:	50 seconds



Item Code: 6300
LB36.0004 Rev B 20160126