



COLOR CHOICES



Red



Blue




FEATURES

- Easy to remove battery pack and compact, lightweight frame
- Conveniently disassembles into four manageable pieces
- Fits in most vehicle trunks
- Top speed of 3.5 mph and enough power to get you where you want to go



The First Name in Travel Mobility®



	Go-Go® ES2	Go-Go® Ultra X	Go-Go Elite Traveller®	Go-Go Elite Traveller® Plus	Go-Go® Sport	Go-Go® LX w/CTS Suspension	Go-Chair®
Max weight capacity:	250 lbs.	260 lbs.	300 lbs.	300 lbs.	325 lbs.	300 lbs.	250 lbs.
Maximum speed: ¹	Up to 4 mph	Up to 4 mph	Up to 4 mph	Up to 4.5 mph	Up to 4 mph	Up to 4.4 mph	Up to 3.5 mph
Ground clearance: ^{2,6}	1.5" at mid-frame	3W: 2" / 4W: 2"	3"	3"	3"	1.5" at motor	1.75"
Turning radius: ²	31"	3W: 33" / 4W: 44"	3W: 34" / 4W: 45.75"	3W: 34" / 4W: 45.75"	3W: 36" / 4W: 52"	3W: 33" / 4W: 44.75"	26.75" without foot platform
Overall length: ²	37"	3W: 37" / 4W: 39.5"	3W: 40" / 4W: 40.25"	3W: 40" / 4W: 40.25"	43"	3W: 39.75" / 4W: 41"	29" without foot platform
Overall width: ²	19.5"	19.5"	21"	21"	20"	20.5"	19"
Battery charger:	Off-board, 1.5 amp	Off-board, 2 amp	Off-board, 2 amp	Off-board, 2 amp	Off-board, 2 amp	Off-board, 2 amp	Off-board, 2 amp
Per-charge range: ^{1,5}	Up to 6 miles	Up to 8 miles	Up to 7.7 miles (w/ 2x 12AH)	Up to 13.3 miles	3W: Up to 9 miles / 4W: Up to 12 miles	3W: Up to 9.4 miles / 4W: Up to 10.5 miles	Up to 8 miles
Battery requirements: ³	(2) 10AH	(2) 12AH	(2) 12AH or (2) 18AH	(2) 18AH	(2) 18AH	(2) 12AH or (2) 18AH	(2) 12 AMP
Battery pack weight: ⁴	1.2 lbs.	12 AH: 18 lbs. 18 AH: 28 lbs.	12 AH: 18 lbs. 18 AH: 28 lbs.	28 lbs.	28 lbs.	18 lbs. / 28 lbs.	20 lbs.
Weight w/o battery pack: ²	58.4 lbs.	3-wheel: 75 lbs. 4-wheel: 78.5 lbs.	3-wheel: 75 lbs. 4-wheel: 86.5 lbs.	3-wheel: 84 lbs. 4-wheel: 86.5 lbs.	3-wheel: 83.2 lbs. 4-wheel: 99.2 lbs.	3-wheel: 85.4 lbs. 4-wheel: 87.4 lbs.	82 lbs.
Heaviest piece: ²	42 lbs. (scooter base)	3-wheel: 29.5 lbs. (front section) 4-wheel: 34 lbs. (front section)	3-wheel: 29.5 lbs. (front section) 4-wheel: 34 lbs. (front section)	3-wheel: 32.5 lbs. (rear section) 4-wheel: 35 lbs. (front section)	3-wheel: 36.2 lbs. (front section) 4-wheel: 52.2 lbs. (front section)	3-wheel: 34 lbs. (rear section) 4-wheel: 36.6 lbs. (front section)	34.5 lbs. (rear section)

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
 2. Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of (+ or -) 3%.
 3. AGM or gel cell type recommended.
 4. Battery weight may vary based on manufacturer.
 5. Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
 6. Measured at midframe.
 7. All shroud panels included.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

The First Name in Travel Mobility®
Pride®
 Mobility Products Corp.
 — LIVE YOUR BEST! —

GOGO[®] es²

COLOR CHOICE



Black



FEATURES

- Two speed settings - high and low (indoor/outdoor)
- Modular design for easy serviceability
- Removable front canvas basket
- Lightweight, one-piece frame for added stability and easy portability

GOGO[®] ULTRA X

ALSO AVAILABLE IN 3-WHEEL.

COLOR CHOICES



Red



Blue



FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Pride's exclusive black, non-scuffing tires
- Auto-connecting front to rear cable
- All-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage

GOGO[®] ELITE TRAVELLER[®]

ALSO AVAILABLE IN 4-WHEEL.

COLOR CHOICES



Red



Blue



FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Pride's exclusive black, non-scuffing tires
- Modular design for easy serviceability
- All-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage

GOGO[®] ELITE TRAVELLER PLUS[®]

ALSO AVAILABLE IN 3-WHEEL.

COLOR CHOICES



Red



Blue



FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Standard 18 AH battery pack
- Front tiller-mounted basket
- All-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage

GOGO[®] SPORT

ALSO AVAILABLE IN 3-WHEEL.

COLOR CHOICES



Red



Blue



FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Pride's exclusive black, non-scuffing tires
- Charger port in tiller
- All-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage

GOGO[®] LX WITH CTS SUSPENSION

ALSO AVAILABLE IN 4-WHEEL.

COLOR CHOICES



Red



Blue

OPTIONAL



Orion Yellow



Onyx Black



Pearl Pink



Sunset Orange



INTERCHANGEABLE SHROUDS

FEATURES

- CTS Suspension (Comfort-Trac front and rear independent suspension)
- Feather-touch disassembly permits simple frame separation with only one hand
- All-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage