

MAXICOMFORT™ CLOUD SERIES

Rise

Available NOW!



Cirrus Medium PR-508M



Cirrus Small PR-508S



Cirrus Large PR-508L

Softouch
with STAIN DEFENSE




Palomino




Admiral




Evergreen



Cabernet



Pearl



Brisa Buckskin*
* upcharge applies

Featuring the Industry's First Pocket Coil Seat Springs!



MAXI COMFORT™
Ultimate Recline Technology

Cirrus
by Golden

The all-new Cirrus™ with MaxiComfort lift chair has taken comfort to a whole new level! The industry's first lift chair featuring pocket coil seat springs, the Cirrus joins the Cloud PR510 as the most comfortable chairs we've ever built! Our pocket coil spring technology and supportive foam casing molds to your lower body, gently cradling you in plush comfort. Experience the newest heights in lift chair comfort – experience the Cirrus!




Rise. Roll. Ride. GOLDEN technologies
YOUR LIFE IN MOTION

Cirrus™



MAXICOMFORT™ CLOUD SERIES

| Model Number | PR-508S | PR-508M | PR-508L |
|--|--------------|--------------|--------------|
| Size | Small | Medium | Large |
| Weight Capacity | 300 LB | 375 LB | 375 LB |
| Recline Positions | Maxi Comfort | Maxi Comfort | Maxi Comfort |
| Standard Fabrics ¹ | 6 | 6 | 6 |
| Fabric Options | All | All | All |
| Chaise Seat | Yes | Yes | Yes |
| Back Type | Biscuit | Biscuit | Biscuit |
| Head Rest Cover | Yes | Yes | Yes |
| Arm Covers | No | No | No |
| Battery Back Up | Yes | Yes | Yes |
| Overall Width | 31.5" | 32" | 35" |
| Overall Height | 41" | 43" | 46" |
| Width Between Arms | 19.5" | 20" | 22" |
| Floor to Top of Seat | 19.5" | 20.5" | 21.5" |
| Seat Depth | 18" | 20" | 21" |
| Seat to Top of Back | 27" | 29" | 31" |
| Seat to Top of Arm | 5.5" | 5.5" | 6.5" |
| Seat Angle Raised | 15° | 15° | 16° |
| Floor to Top of Foot Rest Reclined | 23" | 23" | 24" |
| Floor to Top of Back Reclined | 23" | 22" | 23" |
| Back Angle Closed | 100° | 100° | 105° |
| Back Angle Reclined | 170° | 170° | 175° |
| Foot Rest Angle Closed | 90° | 90° | 90° |
| Foot Rest Angle Opened | 180° | 180° | 180° |
| Distance Required from Wall Reclined | 24" | 26" | 28" |
| Extension Reclined from Back of Seat to End of Foot Rest | 40" | 42" | 43" |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism | Yes | Yes | Yes |
| Warranty on Electrical parts | 3 YR | 3 YR | 3 YR |
| Pro Rated Years 4 Thru 7 | Yes | Yes | Yes |
| Warranty on Mechanical Labor | 3 YR | 3 YR | 3 YR |
| Warranty on Electrical Labor | 1 YR | 1 YR | 1 YR |
| Removable Back Option Available | No | No | No |
| Smart Tek System Standard | Yes | Yes | Yes |
| Customizable * | No | No | No |
| Left Hand Control Option Available | Yes | Yes | Yes |
| Extra Pocket Option Available | Yes | Yes | Yes |
| Heat & Massage Option Available | Yes | Yes | Yes |
| Foot Rest Extension Option Available | No | Yes | Yes |
| Straight Lift Option Available | No | No | No |
| User Height (suggested) | 5'0"-5'3" | 5'4"-5'10" | 5'11"-6'2" |

¹ Buckskin Brisa is a standard, stock fabric on PR508M. Upcharge applies on all sizes.

* See Customize Chart for Dimension limits
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.

Hand-blown, Overstuffed, Conjugated Microfiber Dacron Fill
SUPER SUPPORTIVE

Quick Release AutoDrive Hand Control
EASY DISCONNECT

Widest Fabric Selection in Industry
FITS IN ANY HOME DÉCOR

Adjustable Anti-Skid/Scrape Leg Levelers
STABILITY & SUPPORT

Kiln-dried Northern Maple Hardwood
 • Double-doweled at Stress Points
 • Reinforced Miter Blocks
 • Stapled and Glued
STRONGEST IN THE INDUSTRY

10 Gauge Steel Scissor Mechanism
STRONG & LONG LASTING


Rear Landing Legs
EXTRA STABILITY

5 Heavy Duty Springs
LONG LASTING FOR INCREDIBLE COMFORT



PR508M shown

Our patented technologies, combined with our superior kiln-dried northern maple hardwood frame design, provide the most comfortable lift and recline experience.



The Cirrus™ features a combination of pocket coil seat springs surrounded by foam that molds to your lower body for the ultimate support and comfort.



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

Our luxurious lift and recline chairs are equipped with electronics to help assure a secure lift every time you need it.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68
 Listing # E112417
 FNAE 80-214A

