


V A L U E S E R I E S

Rise



Monarch Medium
PR-355M

SoftTouch
with STAIN DEFENSE



Autumn




Cornflower



Spring



Rosewood



Chestnut*

Std. only on
PR200, PR355M,
PR355L



Capri
PR-200



Monarch Large
PR-355L



Monarch Plus Medium
PR-359M



Monarch Plus Large
PR-359L

Value Series

The Value Series offers three well appointed models at a great value that offer hours of comfortable relaxation. The Capri features a classic, button-back design, while the Monarch models have a seam back. The Capri and Monarch have a split chaise pad and the Monarch Plus models feature a full chaise pad.



Rise. Roll. Ride.



YOUR LIFE IN MOTION



	Capri	Monarch	Monarch	Monarch Plus	Monarch Plus
Model Number	PR-200	PR-355M	PR-355L	PR-359M	PR-359L
Size	Medium	Medium	Large	Medium	Large
Weight Capacity	325 Lb.	375 Lb.	375 Lb.	375 Lb.	375 Lb.
Recline Positions	2	3	3	3	3
Standard Fabrics	5	5	5	4	4
Fabric Options	NONE	ALL	ALL	ALL	ALL
Chaise Seat	NO	NO	NO	YES	YES
Back Type	Button	Seam	Seam	Seam	Seam
Head Rest Cover	YES	YES	YES	YES	YES
Arm Covers	YES	YES	YES	YES	YES
Battery Back Up	YES	YES	YES	YES	YES
Overall Width	31"	32"	34"	32"	34"
Overall Height	41"	42"	44"	41.5"	43"
Width Between Arms	20"	21"	23"	21"	23"
Floor to Top of Seat	18"	18.5"	18.5"	19.5"	20.5"
Seat Depth	19"	21.5"	23"	21.5"	22.5"
Seat to Top of Back	27"	28"	30"	27.5"	29.5"
Seat to Top of Arm	7"	7"	7.5"	6"	7.5"
Floor to Back Edge of Seat Raised	27"	27"	27"	28"	29"
Floor to Front Edge of Seat Raised	22"	23"	23"	24"	25"
Seat Angle Raised	16°	16°	16°	16°	16°
Floor to Top of Foot Rest Reclined	23"	23"	22"	23"	23"
Floor to Top of Back Reclined	35"	33"	33"	33"	33"
Back Angle Closed	110°	110°	110°	110°	110°
Back Angle Reclined	135°	140°	140°	140°	140°
Foot Rest Angle Closed	90°	90°	90°	90°	90°
Foot Rest Angle Opened	160°	176°	176°	176°	176°
Distance Required from Wall Reclined	11"	17"	18"	17"	18"
Extension Reclined from Back of Seat to End of Foot Rest	40"	43"	43"	43"	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES	YES	YES
Warranty on Electrical Parts	1 Yr.	2 Yrs.	2 Yrs.	2 Yrs.	2 Yrs.
Pro Rated Years 4 Thru 7	NO	NO	NO	NO	NO
Warranty on Mechanical Labor	6 Mos.	6 Mos.	6 Mos.	6 Mos.	6 Mos.
Warranty on Electrical Labor	6 Mos.	6 Mos.	6 Mos.	6 Mos.	6 Mos.
Removable Back Option Available	NO	YES	YES	YES	YES
Smart Tek System Standard	YES	YES	YES	YES	YES
Customizable	NO	NO	NO	NO	NO
Left Hand Control Option Available	YES	YES	YES	YES	YES
Extra Pocket Option Available	YES	YES	YES	YES	YES
Heat and Massage Option Available	NO	YES	YES	YES	YES
Foot Rest Extension Option Available	NO	NO	NO	YES	YES
Straight Lift Option Available	YES	YES	YES	YES	YES

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.

Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Battery Backup System



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68 Listing # E112417 FNAE 80-214A

