



### Home Bed Assist Rail Folding Bed Board Combo

Item #	UOM
15062	1/cs

- Assistance for getting in and out of bed. Can be used on either side of the bed. (Figure A)
- Tri-fold bed board for easy portability and storage. (Figure B)
- Removable handle.
- Easy, quick, no tool assembly. Comes in two pieces.
- Simply slides under the mattress for extra-firm support. (Figure C)
- Removable powder coated steel handle enables user to remove when not in use.
- Handle with the "Mid Bar", creates grip for any height use.
- Retail Packaged.

specifications	
Rail Dimensions	11.5" (W) x 22.5" (H) x 30" (D)
Weight Limit	250 lbs.
Warranty	Limited Lifetime

### Home Bed Assist Handle

Item #	UOM
15064	1/cs

- Assistance for getting in and out of bed. Can be used on either side of the bed. (Figure A)
- Removable handle.
- Easy, quick, no tool assembly. Comes in two pieces.
- Safety strap wraps around mattress or box spring to ensure a secure fit.
- Special extra-long safety strap enables user to be in a single, twin, double, queen or king size bed.
- Made of powder coated steel for long life and easy cleaning.
- Handle with multiple crossbars creates easy grip for any height use.
- Retail Packaged.

specifications	
Rail Dimensions	11.5" (W) x 20.75" (H) x 21.5" (D)
Weight Limit	250 lbs.
Warranty	Limited Lifetime